

MOLBIO DIAGNOSTICS PRIVATE LIMITED ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

Mfg. Lic. No.: MFG/IVD/2019/000003

Date: 22/02/2021

PRODUCT: Trueprep® AUTO MTB Sample Pre-treatment Pack

Trueprep® AUTO MTB Lot No.: MSPT089 Mfg Date: 2021-01 Exp Date: 2022-12

Sample Pre-treatment Pack A.R. No.: B-MSPT-089

Pack size: 50T, 20T, 5T STORAGE: 2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
Α	Contents in the carton (MTB Sample Pre-treatment Pack)	(Present (P)/N	ot Present (NP))
1	Lysis Buffer	Present	Present
2	Liquefaction Buffer	Present	Present
3	Disposable Transfer Pipette (1ml)	Present	Present

	Performance	on True	elab [®] Real Time micro PCR An	alyzers		
Sr. No.	Controls	The state of the s			st Lot Values)	
1	Specificity Tests	Test	IC	Test	IC	
Negative	Control					
1	Sample 1	ND	Control should be detected	ND	30.0	
2	Sample 2	ND	Control should be detected	ND	29.0	
3	Sample 3	ND	Control should be detected	ND	29.0	
Positive	Control					
1	Sample 1	Test a	nd Control should be detected	27.0	28.14	
2	Sample 2	Test a	nd Control should be detected	27.33	27.33	
3	Sample 3	Test a	nd Control should be detected	27.17	30.33	

^{*}Ct- Cycle threshold

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALE CHOST

QUALITY CONTROL
Name: Dr. Praveen Kumar
Designation: QC Head

Page 1 of 1

^{*}IC- Internal Control

^{*}ND- Not Detected



MOLBIO DIAGNOSTICS PRIVATE LIMITED ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

Mfg.Lic. No.: MFG/IVD/2020/000119

Date: 22/02/2021

PRODUCT: Trueprep® AUTO Universal Sample Pre-treatment Pack

Trueprep® AUTO Universal Sample Pre-treatment Pack Lot No.: USPT143 Mfg Date: 2021-01 Exp Date: 2022-12

A.R. No.: B-USPT-143a

Pack size: 50T, 20T, 5T

STORAGE: 2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
Α	Contents in the carton (Trueprep® AUTO Universal Sample Pre- treatment Pack)(Present (P)/Not Present (NP))	(Present (P)/ Not	Present (NP))
1	Lysis Buffer	Present	Present
2	Disposable Transfer Pipette (1ml)	Present	Present

	Perform	nance on	Truelab® micro PCR Analyzers	3	
Sr. No.	Control	S	pecification (Ct Values)	Test Lot (Values)	Ct
1	Specificity Tests	Test	IC	Test	IC
Negat	tive Control				
1	Neg HBV	ND	Control should be detected	ND	28.29
2	Neg HIV	ND	Control should be detected	ND	29.75
3	Neg COVID-19	ND	Control should be detected	ND	24.33
Positi	ve Control				-
1	Pos HBV	Test ar	nd Control should be detected	28.33	28.0
2	Pos HCV	Test ar	nd Control should be detected	29.7	30.33
3	Pos COVID-19	Test ar	nd Control should be detected	27.75	24.4

^{*}Ct- Cycle threshold

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS GNOST

PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALE

QUALITY CONTRO Name: Dr. Praveen Kumar **Designation:QC Head**

Page 1 of 1

Molbio Diagnostics Private Limited

^{*}IC- Internal Control

^{*}ND- Not Detected



MOLBIO DIAGNOSTICS PRIVATE LIMITED ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

Mfg.Lic. No.: MFG/IVD/2020/000119

Date: 20/02/2021

PRODUCT: Trueprep® AUTO Universal Sample Pre-treatment Pack

Trueprep® AUTO Universal Lot No.: USPT142 Mfg Date: 2021-01 Exp Date: 2022-12

Sample Pre-treatment Pack A.R. No.: B-USPT-142b

Pack size: 50T, 20T, 5T STORAGE: 2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
A	Contents in the carton (Trueprep® AUTO Universal Sample Pre- treatment Pack)(Present (P)/Not Present (NP))	(Present (P)/ Not	Present (NP))
1	Lysis Buffer	Present	Present
2	Disposable Transfer Pipette (1ml)	Present	Present

	Perfor	mance on	Truelab® micro PCR Analyzers	;	
Sr. No.	Control	S	pecification (Ct Values)	Test Lot (Values)	Ct
1	Specificity Tests	Test	IC	Test	IC
Negati	ive Control				
1	Neg HBV	ND	Control should be detected	ND	27.2
2	Neg HIV	ND	Control should be detected	ND	30.75
3	Neg COVID-19	ND	Control should be detected	ND	22.2
Positiv	ve Control				
1	Pos HBV	Test ar	d Control should be detected	29.8	28.67
2	Pos HCV	Test ar	d Control should be detected	28.07	27.2
3	Pos COVID-19	Test ar	d Control should be detected	29.13	21.0

^{*}Ct- Cycle threshold

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALES

Name: Dr. Prayeen Kumar Designation: QC Head

Page 1 of 1

Molbio Diagnostics Private Limited

^{*}IC- Internal Control

^{*}ND- Not Detected



MOLBIO DIAGNOSTICS PRIVATE LIMITED ISO 13485:2016 certified company CERTIFICATE OF ANALYSIS

Mfg.Lic. No.: MFG/IVD/2019/000003

Date: 20/02/2021

Product: Trueprep® AUTO Universal Cartridge Based Sample Prep Kit

Trueprep® AUTO Universal

Pack size: 5, 25, 50 TESTS

Lot No.: USPK442-P Mfg Date: 2020-09 Exp Date: 2022-08

Cartridge Based

A.R. No.: B-USPK-442-P

Sample Prep Kit

Storage: 2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
Α	Contents in the carton (Trueprep® AUTO Universal Cartridge Based Sample prep kit)	(Present (P)/Not	Present (NP))
1	Universal Reagent Pack	Present	Present
2	Universal Cartridge Pack	Present	Present
3	Disposable Transfer Pipette (3 ml)	Present	Present
4	Package Insert	Present	Present

	Performance on Truelab® micro PCR Analyzers				
Sr. No	Sample Negative/Positive	Negative: Test- ND Positive: Test- Detected	IC- Detected IC- Detected		
		Test (Ct)	IC (Ct)		
1	Negative MTB	ND	29.5		
2	Positive MTB	27.33	28.8		
3	Negative HBV	ND	29.0		
4	Positive HBV	30.13	27.2		
5	Negative HIV	ND	28.0		
6	Positive HCV	27.88	30.33		

^{*}Ct- Cycle threshold

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SAMPLES TESTED AND LOT IS RELEASED.

QUALITY CONTROL Name: Dr. Praveen Kumar Designation: QC Head

VERNA

Page 1 of 1

^{*}IC- Internal Control

^{*}ND- Not Detected