


# GLOBAL FIRST - RAPID POINT-OF-CARE MDx SYSTEM



• Real Time Quantitative **micro** PCR System •

**Sample to Results in less than 1 Hour!**

**Trueprep<sup>®</sup>**  
**AUTO v2**




## **True lab<sup>®</sup> Uno Dx Real Time Quantitative micro PCR Analyzer**

Fully automatic Real Time Quantitative micro PCR Analyzer, Three wavelength system, performs between 10-12 tests in 8 hours. Controls amplification, detection and reporting of results.

## **True lab<sup>®</sup> Duo Real Time Quantitative micro PCR Analyzer**

Fully automatic Real Time Quantitative micro PCR Analyzer, Two channel-Three wavelength system, performs 20-24 tests in 8 hours.



## **True lab<sup>®</sup> Quattro Real Time Quantitative micro PCR Analyzer**

Fully automatic Real Time Quantitative micro PCR Analyzer, Four channel-Three wavelength system, performs 40-48 tests in 8 hours.

## **Trueprep<sup>®</sup> AUTO v2 Universal Cartridge based Sample Prep Device**

The fully automatic sample prep device works in tandem with **Trueprep<sup>®</sup> AUTO** cartridge and **Trueprep<sup>®</sup> AUTO v2** Reagent Kit for extraction and purification of nucleic acids from clinical specimen.



## **True lab<sup>®</sup> micro PCR printer**

Bluetooth printer, prints wirelessly the results of the PCR tests performed by the **True lab<sup>®</sup> UnoDx/Duo/Quattro Real Time Quantitative micro PCR Analyzer**

## **Truepet<sup>®</sup> SPA Fixed volume (6µl) Precision micropipette**

Pre-calibrated Single Push Auto ejector fixed volume (6µl) micropipette that ensures accurate pipetting of elute.



# ANY WHERE, ANY TIME, REAL TIME PCR

## Truenat<sup>®</sup> Real Time micro PCR Chips

### PRODUCTS AVAILABLE


- |                                                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                  |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• MTB <sup>CE</sup></li> <li>• MTB Plus <sup>CE</sup></li> <li>• MTB-RIF Dx <sup>CE</sup></li> <li>• MTB-INH</li> <li>• HAV <sup>CE</sup></li> <li>• HBV with VL</li> <li>• HCV with VL</li> <li>• HEV <sup>CE</sup></li> <li>• HIV-1 with EID+VL</li> <li>• HIV-1/HIV-2</li> </ul> | <ul style="list-style-type: none"> <li>• H1N1 <sup>CE</sup></li> <li>• H3N2/H1N1</li> <li>• Chlamydia</li> <li>• Gonorrhoea <sup>CE</sup></li> <li>• Chlamydia/Gonorrhoea</li> <li>• Trichomonas <sup>CE</sup></li> <li>• Dengue <sup>CE</sup></li> <li>• Chikungunya <sup>CE</sup></li> <li>• Dengue/Chikungunya <sup>CE</sup></li> <li>• Malaria Pf <sup>CE</sup></li> </ul> | <ul style="list-style-type: none"> <li>• Malaria Pv/Pf <sup>CE</sup></li> <li>• Salmonella <sup>CE</sup></li> <li>• Leptospira <sup>CE</sup></li> <li>• Shigella</li> <li>• Scrub Typhus <sup>CE</sup></li> <li>• Rabies <sup>CE</sup></li> <li>• HLA-B27 <sup>CE</sup></li> <li>• HPV-HR <sup>CE</sup></li> <li>• Nipah <sup>CE</sup></li> <li>• GBS <sup>CE</sup></li> </ul> | <ul style="list-style-type: none"> <li>• Influenza A/B <sup>CE</sup></li> <li>• Beta CoV <sup>CE</sup></li> <li>• COVID-19 <sup>CE</sup></li> <li>• SARS CoV-2 <sup>CE</sup></li> <li>• Clostridium difficile</li> <li>• HSV (1/2)</li> <li>• Cholera</li> </ul> |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

### PRODUCTS COMING SOON

- |                                                                                                                                                                                                                                                                                            |                                                                                                                                                                                                                                                                                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• HPV-HR 16 &amp; 18</li> <li>• HPV-HR 16,18/31,39,45</li> <li>• MTB Ultima</li> <li>• Staph/MRSA</li> <li>• ZIKA</li> <li>• MTB-FQ</li> <li>• MTB/NTM</li> <li>• CMV</li> <li>• Mumps</li> <li>• Dengue/Zika</li> <li>• Rubella/Measles</li> </ul> | <ul style="list-style-type: none"> <li>• Staph/PA</li> <li>• Influenza AB/COVID-19</li> <li>• Rotavirus</li> <li>• Mucormycosis</li> <li>• Mycobacterium leprae</li> <li>• Monkeypox (MPX)</li> <li>• Leishmania</li> <li>• Mycoplasma genitalium</li> <li>• MTB Ultima/COVID-19</li> <li>• EBV</li> <li>• Brucella</li> </ul> | <ul style="list-style-type: none"> <li>• HHV-6 (A/B)</li> <li>• Syphilis</li> <li>• Pseudomonas aeruginosa (PA)</li> <li>• Streptococcus pneumoniae (SP)</li> <li>• Kysanur Forest Disease (KFD)</li> <li>• Haemophilus influenzae B (HIB)</li> <li>• Neisseria meningitidis (Nm)</li> <li>• Japanese encephalitis virus (JEV)</li> <li>• MERS (Middle East Respiratory Syndrome)</li> <li>• Respiratory syncytial virus (RSV)</li> <li>• Sickle Cell Anemia screening</li> </ul> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|


### PRODUCTS PIPELINE

- |                                                                                                                                                |                                                                                                                                                                                          |                                                                                                                                      |
|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• HDV</li> <li>• Ebola</li> <li>• Adenovirus</li> <li>• Salmonella AMR</li> <li>• Klebsiella</li> </ul> | <ul style="list-style-type: none"> <li>• Bordetella pertussis</li> <li>• H5N1</li> <li>• Rift valley fever</li> <li>• Tick Borne encephalitis</li> <li>• Rickettsial diseases</li> </ul> | <ul style="list-style-type: none"> <li>• Yellow fever virus</li> <li>• Listeria monocytogens</li> <li>• Toxoplasma gondii</li> </ul> |
|------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|



### Truenat<sup>®</sup> Chips

Disease specific, ready-to-use and disposable micro PCR chips that run on the **Truelab<sup>®</sup> Uno Dx/ Duo/Quattro** to get quantitative real time PCR results.



### Trueprep<sup>®</sup> AUTO Cartridge

Disposable plastic cartridge for automated extraction of nucleic acids from clinical specimens.

### Trueprep<sup>®</sup> AUTO v2 Reagent Kit

Reagent kit for use with **Trueprep<sup>®</sup> AUTO v2** cartridge for automated processing of clinical specimens



### Truelab<sup>®</sup> Real Time PCR Workstation Field Case

Compact Portable carry case.

\*Available only with **Truelab<sup>®</sup> Uno Dx**.



**Molbio Diagnostics Private Limited**

**Manufacturing & Sales Office:** Plot No. L-46, Phase II D, Verna Industrial Estate, Verna, Goa-403722, INDIA

Tel: +91-832-2783267 **Email:** sales@molbiodiagnostics.com, customersupport@molbiodiagnostics.com

**Website:** www.molbiodiagnostics.com

**ANY WHERE, ANY TIME, REAL TIME PCR**