

**MOLBIO DIAGNOSTICS LIMITED**  
**EN ISO 13485:2016 certified company**  
**CERTIFICATE OF ANALYSIS**

**License No.:** MFG/IVD/2019/000003

**Date:** 10/06/2025

**Doc ID:** COA/TP/ATM106/50T

<b>Product Name</b>	Trueprep <sup>®</sup> AUTO Transport Medium for Swab Specimen Pack
<b>Generic Name</b>	Medium for Collection, Decontamination and Transport of Swab Specimen for Trueprep <sup>®</sup> AUTO
<b>Lot No.</b>	ATM106
<b>Mfg. Date</b>	2025-05
<b>Exp. Date</b>	2027-04
<b>Pack Size</b>	50T
<b>Storage condition</b>	2°C to 40°C

**PRECAUTION:** Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
<b>A</b>	<b>Contents in the carton (Trueprep<sup>®</sup> AUTO Transport Medium for Swab Specimen Pack )</b>	<b>(Present (P)/Not Present (NP))</b>	
1	Transport Medium for Swab Specimen Tubes	Present	Present
2	Package Insert	Present	Present

Performance on Truelab® Real Time micro PCR Analyzer						
Sr. No.	Control		Specification (Ct Values)		Test Lot (Ct Values)	
	Specificity Tests	Test	IC	Test 1 E.gene (Average)	Test 2 Orf1A (Average)	IC (Average)
Negative Control						
1	Negative COVID-19	ND	Control should be detected	ND	ND	27.74
Positive Control						
1	Positive COVID-19	Test and Control should be detected		30.31	31.06	26.7
Non Template Control Test for COVID-19						
1	NTC	Test 1, Test 2 and Control should not be detected		ND	ND	ND

\*Ct-Cycle Threshold \*IC- Internal Control \*ND- Not Detected \*NA- Not Applicable \*NTC- Non Template Control

The lot is tested and verified as per QC Standards and is approved.

   
Name: Dr. Praveen Kumar MK  
Department: Quality Control  
Designation: QC Head