

MOLBIO DIAGNOSTICS LIMITED

EN ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

Mfg. Lic. No.: MFG/IVD/2019/000003

Date: 13/11/2025

Doc ID: COA/TN/PCIV017

BRAND NAME	Truenat [®] Positive Control Kit - Panel IV
GENERIC NAME	Control Kit for Truenat [®] Chip-based micro PCR Tests: HAV, HEV, Rabies, LTS, Nipah and Malaria Pv/Pf
LOT NO	PCIV017
MFG. DATE	2025-10
EXP. DATE	2026-09
PACK SIZE	08T
STORAGE	2°C to 40°C A.R. No. : B-PCIV-017

PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
1	Negative control	Present	Present
2	Positive control (08- strip tubes containing dried down positive control)	Present	Present
3	Package insert	Present	Present

Observation (Ct Value)			
Sr. No.	Truenat [®] test name	Positive control	Negative control
1	HAV	27.69 ± 0.51	ND
2	HEV	27.07 ± 0.55	ND
3	Rabies	27.65 ± 0.84	ND
4	LTS	27.30 ± 0.26	ND
5	Nipah	27.44 ± 0.49	ND
6	Malaria Pv/Pf	Pv: 26.68 ± 0.72 Pf: 25.62 ± 0.54	Pv: ND Pf: ND

*ND- Not detected, Ct- Cycle threshold

*A deviation of =<1.5 is acceptable for Ct values

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALE

fa *km*
QUALITY CONTROL ★
Name: Dr. Praveen Kumar MK
Designation: QC Head

Page 1 of 1