

**MOLBIO DIAGNOSTICS PRIVATE LIMITED**  
**ISO 13485:2016 certified company**  
**CERTIFICATE OF ANALYSIS**

Mfg. Lic. No.: MFG/IVD/2019/000003

Date: 12/12/2025  
Doc ID: COA/TN/PCV002

BRAND NAME	Truenat <sup>®</sup> Positive Control Kit - Panel V
GENERIC NAME	Control Kit for Truenat <sup>®</sup> Chip-based micro PCR Tests: H3N2, Zika, Leish, Scrub T, Cholera and Shigella
LOT NO	PCV007
MFG. DATE	2025-12
EXP. DATE	2026-11
PACK SIZE	08T
STORAGE	2°C to 40°C
A.R. No. : B-PCV-007	

**PRECAUTION:** Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
1	Negative control	Present	Present
2	Positive control (08- strip tubes containing dried down positive control)	Present	Present
3	Package insert	Present	Present

Observation (Ct Value)			
Sr. No.	Truenat <sup>®</sup> test name	Positive control	Negative control
1	H3N2	25.20 ± 0.35	ND
2	Zika	27.35 ± 0.32	ND
3	Leish	27.83 ± 0.21	ND
4	Scrub T	27.08 ± 0.42	ND
5	Cholera	27.52 ± 0.44	ND
6	Shigella	29.89 ± 0.51	ND

\*ND- Not detected, Ct- Cycle threshold

**MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCT IS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALE**

for   
QUALITY CONTROL  
Name: Dr. Praveen Kumar MK  
Designation: QC Head



Page 1 of 1