

MOLBIO DIAGNOSTICS LIMITED
EN ISO 13485:2016 certified company
CERTIFICATE OF ANALYSIS

Licence No.: MFG/IVD/2019/000018

Date: 21/02/2026

Document ID:COA/I000TN/I0001TX

| | |
|-------------------|--|
| Product Name | Truenat [®] MTB-RIF Dx |
| Generic Name | Chip-based Real Time PCR Test for Rifampicin Resistant <i>Mycobacterium tuberculosis</i> |
| Lot No. | I0001TX |
| Set No. | NA |
| Mfg. Date | 2026-01 |
| Exp. Date | 2027-12 |
| Pack Size | 5T, 20T, 50T |
| Storage condition | 2°C to 30°C |

PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

| Sr. No. | Components in each kit | Specification | Observation |
|---------|--|--|-------------------------|
| 01 | Number of pouches as per pack size | Should be present | Present |
| 02 | Pouches Contents: 1 chip, 1 micro tube, 1 desiccant pouch and 1 filter barrier tip | Should be present | Present |
| 03 | Colour of desiccant pouch | Should be present and Orange in Colour | Present / Orange Colour |
| 04 | Package Insert | Should be present | Present |
| 05 | Seal integrity test | Should pass | Pass |

| Performance testing data | | | | |
|--------------------------|--------------|--------------|--------------|-----------------------------|
| Samples | Ct Values | | | Results |
| | Tm | | Control | |
| | P2 | P1 | | |
| Positive sample 1 | 66.65 ± 0.34 | 67.02 ± 0.45 | 70.25 ± 0.31 | RIF Resistance detected |
| Positive sample 2 | 72.24 ± 0.41 | 66.86 ± 0.48 | 70.39 ± 0.33 | RIF Resistance not detected |
| Positive sample 3 | 72.22 ± 0.34 | 66.97 ± 0.34 | 70.57 ± 0.40 | RIF Resistance not detected |

Tm- Melting temperature, Ct- Cycle threshold

*Mfg- Manufacturing, Exp- Expiry

The lot is tested and verified as per QC Standards and is approved.



Name : Dr. Praveen Kumar MK
Department: Quality Control
Designation: QC Head

for