

**MOLBIO DIAGNOSTICS LIMITED**  
**EN ISO 13485:2016 certified company**  
**CERTIFICATE OF ANALYSIS**

Licence no.: MFG/IVD/2025/000148

Date: 20/03/2026  
Document ID: COA/TN/PCI063

<b>Product Name</b>	Truenat <sup>®</sup> Positive Control Kit - Panel I
<b>Generic Name</b>	Control Kit for Truenat <sup>®</sup> Chip-based micro PCR Tests: MTB, MTB Plus, MTB-RIF Dx, Beta CoV, SARS CoV-2, COVID-19 and H1N1
<b>Lot No.</b>	PCI063
<b>Mfg. Date</b>	2026-03
<b>Exp. Date</b>	2027-02
<b>Pack Size</b>	08T
<b>Storage condition</b>	2°C to 40°C

**PRECAUTION:** Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
01	Negative control	Should be present	Present
02	Positive control (08- strip tubes containing dried down positive control)	Should be present	Present
03	Package insert	Should be present	Present

**MOLBIO DIAGNOSTICS LIMITED**  
**EN ISO 13485:2016 certified company**  
**CERTIFICATE OF ANALYSIS**

Licence no.: MFG/IVD/2025/000148

Date: 20/03/2026  
Document ID: COA/TN/PCI063

Performance testing data			
Sr. No.	Truenat <sup>®</sup> test name	Positive control	Negative control
1	MTB	29.43 ± 0.81	ND
2	MTB Plus	28.97 ± 0.33	ND
3	MTB-RIF Dx	RIF Resistance detected	Indeterminate
4	SARS CoV-2	28.47 ± 0.18	ND
5	COVID-19	Orf 1A: 28.01 ± 0.37 E gene: 27.87 ± 0.49	ND
6	H1N1	28.12 ± 0.90	ND

\*Mfg- Manufacturing, Exp- Expiry

**The lot is tested and verified as per QC standards and is approved.**

  
  
 for Name: Dr. Praveen Kumar MK  
 Department: Quality Control  
 Designation: QC Head