

MOLBIO DIAGNOSTICS LIMITED
EN ISO 13485:2016 certified company
CERTIFICATE OF ANALYSIS

Licence No.: MFG/IVD/2019/000018

Date: 15/04/2026

Document ID: COA/TN/I0003VD

| | |
|-------------------|--|
| Product Name | Truenat [®] SARS CoV-2 |
| Generic Name | Chip-based Real Time PCR Test for COVID-19 |
| Lot No. | I0003VD |
| Set No. | NA |
| Mfg. Date | 2026-03 |
| Exp. Date | 2028-02 |
| Pack Size | 5T, 25T |
| Storage condition | 2°C to 30°C |

PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

| Sr. No. | Components in each kit | Specification | Observation |
|---------|--|--|-------------------------|
| 01 | Number of pouches as per pack size | Should be present | Present |
| 02 | Pouches Contents: 1 chip, 1 micro tube, 1 desiccant pouch and 1 DNase and RNase free filter barrier microtip | Should be present | Present |
| 03 | Colour of desiccant pouch | Should be present and Orange in Colour | Present / Orange Colour |
| 04 | Package Insert | Should be present | Present |
| 05 | Seal integrity test | Should pass | Pass |

| Performance testing data | | |
|--------------------------|--------------|--------------|
| Samples | Ct Values | |
| | RdRP | Control |
| Non Template Control | ND | ND |
| Negative Sample | ND | 27.06 ± 0.39 |
| Positive Sample 1 | 22.32 ± 0.17 | 26.90 ± 0.94 |
| Positive Sample 2 | 25.86 ± 0.27 | 27.32 ± 0.51 |
| Positive Sample 3 | 30.71 ± 0.85 | 27.38 ± 0.83 |

ND – Not detected, Ct- Cycle threshold

Mfg- Manufacturing, Exp- Expiry

The lot is tested and verified as per QC Standards and is approved.

for Name : Dr. Praveen Kumar MK
Department : Quality Control
Designation: QC Head

