

MOLBIO DIAGNOSTICS LIMITED

EN ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

Licence no.: MFG/IVD/2025/000148

Date: 17/04/2026
Document ID: COA/TN/PCIV018

Product Name	Truenat [®] Positive Control Kit - Panel IV
Generic Name	Control Kit for Truenat [®] Chip-based micro PCR Tests: HAV, HEV, Rabies, LTS, Nipah and Malaria Pv/Pf
Lot No.	PCIV018
Mfg. Date	2026-04
Exp. Date	2027-03
Pack Size	8T
Storage condition	2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
01	Negative control	Should be present	Present
02	Positive control (08- strip tubes containing dried down positive control)	Should be present	Present
03	Package insert	Should be present	Present

Performance testing data			
Sr. No.	Truenat [®] test name	Positive control	Negative control
1	HAV	28.59 ± 0.43	ND
2	HEV	28.71 ± 0.51	ND
3	Rabies	25.20 ± 0.53	ND
4	LTS	27.85 ± 0.63	ND
5	Malaria Pv/Pf	Pv: 27.25 ± 0.43 Pf: 28.03 ± 0.52	ND
6	Nipah	27.51 ± 0.62	ND

*Mfg- Manufacturing, Exp- Expiry

The lot is tested and verified as per QC standards and is approved.



for Name: Dr. Praveen Kumar MK
Department: Quality Control
Designation: QC Head