

**MOLBIO DIAGNOSTICS LIMITED**  
**EN ISO 13485:2016 certified company**  
**CERTIFICATE OF ANALYSIS**

**License No.:** MFG/IVD/2025/000148

**Date:** 14/04/2026

**Doc ID:** COA/TP/USPT295

<b>Product Name</b>	Trueprep <sup>®</sup> AUTO Universal Sample Pre-treatment Pack
<b>Generic Name</b>	Universal Sample Pre-treatment pack for Non MTB samples
<b>Lot No.</b>	USPT295
<b>Mfg. Date</b>	2026-03
<b>Exp. Date</b>	2028-02
<b>Pack Size</b>	25T, 5T
<b>Storage condition</b>	2°C to 40°C

**PRECAUTION:** Refer Safety Data Sheet and follow instructions for use.


<b>Sr. No.</b>	<b>Components in each kit</b>	<b>Specification</b>	<b>Observation</b>
<b>A</b>	<b>Contents in the carton (Trueprep<sup>®</sup> AUTO Universal Sample Pre-treatment Pack)</b>	<b>[Present (P)/Absent (A)]</b>	
1	Lysis Buffer	Present	Present
2	Disposable Transfer Pipette (1 ml Graduated)	Present	Present
3	Package Insert	Present	Present

Performance on Truelab <sup>®</sup> Real Time micro PCR Analyzers						
Sr. No.	Specificity Tests	Specification		Test Lot		
		Test	Control	Test 1 (Average)	Test 2 (Average)	Control (Average)
<b>Negative Control</b>						
1	Negative HBV	ND	Control should be detected	ND	NA	26.79
2	Negative HCV	ND	Control should be detected	ND	NA	26.32
<b>Positive Control</b>						
1	Positive HBV	Test and Control should be detected		31.14	NA	27.65
2	Positive HCV	Test and Control should be detected		30.70	NA	26.50
<b>Non Template Control Test for COVID-19</b>						
1	NTC	Test 1, Test 2 and Control should not be detected		ND	ND	ND

\*Ct-Cycle Threshold    \*ND- Not Detected    \*NA- Not Applicable

The lot is tested and verified as per QC Standards and is approved.



*for*  Name: Dr. Praveen Kumar MK  
Department: Quality Control  
Designation: QC Head