

MOLBIO DIAGNOSTICS LIMITED
EN ISO 13485:2016 certified company
CERTIFICATE OF ANALYSIS

Licence no.: MFG/IVD/2025/000148

Date: 21/05/2026
Document ID: COA/TN/PCIII023

Product Name	Truenat [®] Positive Control Kit - Panel III
Generic Name	Control Kit for Truenat [®] Chip-based micro PCR Tests: CT, NG, Trich, HPV-HR, GBS and Salmonella
Lot No.	PCIII023
Mfg. Date	2026-05
Exp. Date	2027-04
Pack Size	8T
Storage condition	2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
01	Negative control	Should be present	Present
02	Positive control (08- strip tubes containing dried down positive control)	Should be present	Present
03	Package insert	Should be present	Present

Performance testing data			
Sr. No.	Truenat [®] test name	Positive control	Negative control
1	CT	27.22 ± 0.63	ND
2	NG	26.20 ± 0.62	ND
3	Trich	27.16 ± 0.67	ND
4	HPV _{-HR}	HPV 16/31: 26.54 ± 0.28 HPV 18/45: 27.13 ± 0.16	ND
5	GBS	28.14 ± 0.87	ND
6	Salmonella	28.01 ± 0.19	ND

*Mfg- Manufacturing, Exp- Expiry

The lot is tested and verified as per QC standards and is approved.



for Name: Dr. Praveen Kumar MK
Department: Quality Control
Designation: QC Head