

MOLBIO DIAGNOSTICS LIMITED
EN ISO 13485:2016 certified company
CERTIFICATE OF ANALYSIS

License No.: MFG/IVD/2019/000003

Date: 13/05/2026

Doc ID: COA/TP/ALR007

Product Name	Trueprep [®] AUTO Liquefaction Reagent
Generic Name	Reagent Pack for Liquefaction of samples
Lot No.	ALR007
Mfg. Date	2026-04
Exp. Date	2028-03
Pack Size	50T
Storage condition	2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
A	Contents in the carton (Trueprep[®] AUTO Sputum Liquefaction Reagent)	[Present (P)/Absent (A)]	
01	Liquefaction Reagent	Present	Present

Performance on Truelab® Real Time micro PCR Analyzer					
Sr. No.	Specificity Tests	Specification		Test Lot	
		Test	Control	Test (Average)	Control (Average)
Negative Control					
01	Negative MTB	ND	Control should be detected	ND	27.11
02	Negative MTB	ND		ND	26.87
03	Negative MTB	ND		ND	26.47
Non Template Control					
01	NTC	Test and Control should not be detected		ND	ND

*Ct- Cycle Threshold

*ND- Not Detected

*NTC- Non Template Control

The lot is tested and verified as per QC Standards and is approved.




Name: Dr. Praveen Kumar MK
Department: Quality Control
Designation: QC Head