

MOLBIO DIAGNOSTICS LIMITED

EN ISO 13485:2016 certified company

CERTIFICATE OF ANALYSIS

License No.: MFG/IVD/2025/000148

Date: 08/09/2026

Doc ID: COA/TP/ATM119

Product Name	Trueprep [®] AUTO Transport Medium for Swab Specimen Pack
Generic Name	Medium for Collection, Decontamination and Transport of Swab Specimen for Trueprep [®] AUTO
Lot No.	ATM119
Mfg. Date	2026-08
Exp. Date	2028-07
Pack Size	20T, 50T
Storage condition	2°C to 40°C

PRECAUTION: Refer Safety Data Sheet and follow instructions for use.

Sr. No.	Components in each kit	Specification	Observation
A	Contents in the carton (Trueprep[®] AUTO Transport Medium for Swab Specimen Pack)	[Present (P)/Absent (A)]	
1	Transport Medium for Swab Specimen Tubes	Present	Present
2	Package Insert	Present	Present
B	Leak Test	Specification	Observation
1	Leak Test	Should Pass	Pass


TEL: 02-6724888
www.molbiolab.com

Performance on Truelab® Real Time micro PCR Analyzers						
Sr. No.	Specificity Tests	Specification		Test Lot		
		Target	Control	Target 1 E.gene (Average)	Target 2 Orf1A (Average)	Control (Average)
Negative Control						
1	Negative COVID-19	ND	Control should be detected	ND	ND	24.69
Positive Control						
1	Positive COVID-19	Target and Control should be detected		22.44	22.50	26.00
Non Template Control Test for COVID-19						
1	NTC	Target 1, Target 2 and Control should not be detected		ND	ND	ND

*ND- Not Detected

*NTC- Non Template Control

The lot is tested and verified as per QC Standards and is approved.



for

Name: Dr. Praveen Kumar MK

Department: Quality Control

Designation: QC Head