

MOLBIO DIAGNOSTICS LIMITED
EN ISO 13485:2016 certified company
CERTIFICATE OF ANALYSIS

Licence no.: MFG/IVD/2025/000148

Date: 26/09/2026
Document ID: COA/TN/PCII030

Product Name	Truenat [®] Positive Control Kit - Panel II
Generic Name	Control Kit for Truenat [®] Chip-based micro PCR Tests: HBV, HCV, HIV-1, HLA-B27, Dengue, Chikungunya and Influenza A/B
Lot No.	PCII030
Mfg. Date	2026-09
Exp. Date	2027-08
Pack Size	8T
Storage condition	2°C to 40°C


PRECAUTION: Refer Safety Data Sheet and Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
01	Negative control	Should be present	Present
02	Positive control (08- strip tubes containing dried down positive control)	Should be present	Present
03	Package insert	Should be present	Present

Performance testing data			
Sr. No.	Truenat® test name	Positive control	Negative control
1	HBV	28.20 ± 0.30	ND
2	HCV	26.17 ± 0.20	ND
3	HIV-1	27.63 ± 0.38	ND
4	HLA-B27	26.96 ± 0.72	ND
5	Dengue/Chikungunya	Dengue: 26.82 ± 0.37 Chikungunya: 27.30 ± 0.30	ND
6	Influenza A/B	Influenza A: 28.82 ± 0.58 Influenza B: 29.51 ± 0.75	ND

*Mfg- Manufacturing, Exp- Expiry

The lot is tested and verified as per QC standards and is approved.

for  

Name: Dr. Praveen Kumar M.K.
Department: Quality Control
Designation: QC Head