

Announcement:

Product Name: Truenat™ HLA-B27

Chip-based real time PCR test for HLA-B27

Molbio Diagnostics Pvt. Ltd., announced the launch of new Chip based real time PCR test **Truenat™ HLA-B27**. This test is designed for qualitative detection and diagnosis of human leukocyte antigen (HLA) serotype B27 related to an auto-immune or immune-mediated disease such as juvenile rheumatoid arthritis or ankylosing spondylitis.

Truenat™ HLA-B27
Chip-based Real Time PCR Test for HLA-B27

**DIAGNOSIS AND MANAGEMENT OF
ANKYLOSING SPONDYLITIS - A SILENTLY PROGRESSING DISEASE**

For Product Enquiry Call Us at +91-832-6724888 or visit us at www.molbiodiagnostics.com