


Announcement:

New Launch: Truenat™ HCV

Chip-based real time PCR Test for Hepatitis C Virus

Molbio Diagnostics Pvt. Ltd., announced the launch of new Chip based real time PCR test **Truenat™ HCV**. This RT-PCR test is designed for quantitative detection of Hepatitis C Virus (HCV) RNA in human plasma, serum and whole blood samples. **Truenat™ HCV** aids in the diagnosis and confirmation of HCV infection and in estimation of viral load.



Hepatitis C is Often Hidden.

Hepatitis C

**TIME TO
TEST → TREAT → CURE**

Truenat™ HCV
Chip-based Real Time PCR Test for
Hepatitis C Virus

Innovation → *Impact*

For Product Enquiry Call Us at +91-832-6724888 or for more details visit us at www.molbiodiagnostics.com