



MOLBIO DIAGNOSTICS PVT. LTD.
An ISO 9001 & ISO 13485 certified company

CERTIFICATE OF ANALYSIS

Mfg. Lic. No. : 874/L

Date: 22/11/2018

PRODUCT : Truenat™ MTB-RIF Dx
LOT. NO. : TX021
MFG. DATE : 2018- 11
EXP. DATE : 2020- 10
PACK SIZE : 5T, 20T
STORAGE : 2°C to 30°C

A.R. NO. : B-TX-021a

PRECAUTION: The product is offered without warranty for any other specific use than that specified in the Instructions for Use.

Sr. No.	Components in each kit	Specification	Observation
1	Pouch contains 1 chip, 1 micro tube, 1 desiccant pouch and 1 filter barrier tip	5 nos, 20 nos.	Pass
2	Colour of desiccant pouch	Blue	Pass
3	Package insert	01 no.	Pass

Performance				
Specificity Tests	Tm		Control	Result
	P1	P2		
Sample 1	66.24	71.82	71.63	RIF RESISTANCE NOT DETECTED
	66.71	72.04	71.14	RIF RESISTANCE NOT DETECTED
	66.58	72.05	70.73	RIF RESISTANCE NOT DETECTED
Sample 2	66.06	64.76	70.84	RIF RESISTANCE DETECTED
	66.52	65.57	70.22	RIF RESISTANCE DETECTED

Performance						
Specificity Tests	Specification					
	P1 Tm			P2 Tm		
Positive sample	RIF resistance detected	<65°C>	RIF resistance not detected	RIF resistance detected	<70.5°C>	RIF resistance not detected

MANUFACTURING RECORDS ARE VERIFIED, CHECKED AND PRODUCTS MANUFACTURED AS PER THE MASTER FORMULA, SAMPLES TESTED AND LOT IS RELEASED FOR SALE

QUALITY CONTROL
Name: Dr. Praveen Kumar
Designation: QC Head

Page 1 of 1

Molbio Diagnostics Pvt. Ltd.

Regd. Office: H. No. 13, Sagar Society, Dona Paula, Panaji, N. Goa 403 004. INDIA
Works: Plot Nos. M 46-47, Phase III - B, Verna Industrial Estate, Verna - Goa 403 722. INDIA
Tel: 91 832 2783267/6624783 E mail: admin@molbiodiagnostics.com
Email: sales@molbiodiagnostics.com Website: www.molbiodiagnostics.com
CIN: U33125GA2000PTC002909