

Announcement:

New Launch: Truenat™ LTS

Chip-based Real Time PCR Test for Leptospira

Molbio Diagnostics Private Limited, announces the launch of new Chip-based real time PCR test **Truenat™ LTS**. It is designed for the semi quantitative detection of Leptospira and aids in the diagnosis of Leptospirosis in human blood/serum/plasma, cerebrospinal fluid (CSF) and urine specimen.

Leptospirosis

Early and accurate diagnosis followed by treatment of Leptospirosis can save your life

WARNING!
LEPTOSPIROSIS
HEALTH HAZARD

FRESH WATER STREAMS AND MUD
POSSIBLY POLLUTED WITH BACTERIA

EXERCISE CAUTION

Truenat™ LTS
Chip-based Real Time PCR Test for Leptospira

For Product Enquiry Call Us at +91-832-6724888 or for more details visit us at www.molbiodiagnostics.com