


**NEW
LAUNCH**


Truenat[®] MTB-INH

Chip-based Real Time PCR test for
Isoniazid Resistant
Mycobacterium tuberculosis

- Truenat[®] is the world's first point-of-care molecular diagnostics platform approved by **WHO** and **ICMR** for detecting TB and Rifampicin resistance
- **Molbio Diagnostics** has now Launched Truenat[®] MTB-INH to detect Isoniazid resistance
- Cases of Isoniazid (INH) mono-resistance have **grown significantly**
- Rifampicin susceptible Isoniazid resistant TB (Hr-TB) is the major cause of treatment failure, relapse and further drug resistance in TB patients
- Truenat[®] MTB-INH is the only point-of-care Real-Time PCR test that aids in diagnosing (Hr-TB) as the frontline test



ENSURES TIMELY AND
ACCURATE DIAGNOSIS.
ALLOWS START OF
EFFECTIVE TREATMENT
THE SAME DAY.



FOR PRODUCT ENQUIRY
+91-832-6724888



VISIT US AT
www.molbiodiagnostics.com